
Ein
Fallbuch

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 **ZimVie**

ZimVie DENTAL SOLUTIONS

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
Kim SJ et al. Effect of bone block graft with rh-BMP2 on vertical bone augmentation. Int J Oral Maxillofac Surg (2010) 39:883-888

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
Keith JD et al. Clinical and histologic evaluation of a mineralized block allograft: Results from the developmental period (2001–2004). Int J Periodont Rest (2006) 26:321-327.

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Literaturangaben

1. Leonetti JA, Koup R. Localized maxillary ridge augmentation with a block allograft for dental implant placement: case reports. *Implant Dent* (2003) 12:217-26.
2. Keith JD. Localized ridge augmentation with a block allograft followed by secondary implant placement: A case report. *Int. J Periodont Rest Dent* (2004) 24:11-17.
3. Minichetti JC, D'amore JC, Hong AY, Cleveland DB. Human histologic analysis of mineralized bone allograft (Puros) placement before implant surgery. *J Oral Implantol* (2004) 30:74-82.
4. Petrunaro PS, Amar S. Localized ridge augmentation with allogenic block grafts prior to implant placement: case reports and histologic evaluations. *Implant Dent* (2005) 14:139-48.
5. Schau I, Plöger M, Von Der Haar A. Blockaugmentation ohne Knochenentnahme. *Oralchirurgie Journal* (2005) 4:34-36.
6. Keith JD, Petrunaro P, Leonetti JA, Elwell CW, Zeren KJ, Caputo C, Nikitakis NG, Schopf C, Warner MM. Clinical and histologic evaluation of a mineralized block allograft: Results from the developmental period (2001-2004). *Int J Periodont Rest Dent* (2006) 26:321-327.
7. Plöger M, Schau I. Rekonstruktion von Kieferkammdefekten mittels allogener Spongiosa. *Implantologie Zeitung* (2006) 10:16-18.
8. Giray B, Meral G. Allogenic grafts in oral surgery: clinical findings and follow-up. *Journal of Hacettepe Faculty of Dentistry* (2007) 31:31-37.
9. Morelli T, Neiva R, Wang HL. Human histology of allogeneic block grafts for alveolar ridge augmentation: case report. *Int J Periodont. Rest Dent* (2009) 29:649-56.
10. Plöger M, Schau I. Kieferkammaugmentation mit allogenen Knochenblöcken. *Journal for Continuing Dental Education* (2009) 12:242-249.
11. Kim SG, Park JS, Lim SC. Placement of implant after bone graft using J block allograft. *Implant Dent* (2010) 19:21-8.
12. Peleg M, Sawatari Y, Marx RN, Santoro J, Cohen J, Bejarano P, Malinin T. Use of corticocancellous allogeneic bone blocks for augmentation of alveolar bone defects. *Int J Oral Maxillofac Implants* (2010) 25:153-62.
13. Plöger M, Schau I. Allogener Knochenblock zur präimplantologischen Augmentation. *Journal for Continuing Dental Education* (2010) 13:84-88.
14. Kim SJ, Shin HS, Shin SW. Effect of bone block graft with rhBMP-2 on vertical bone augmentation. *Int J Oral Maxillofac Surg* (2010) 39:883-8.
15. Jacotti M, Wang HL, Fu JH, Zamboni G, Bernardello F. Ridge augmentation with mineralized block allografts: clinical and histological evaluation of 8 cases treated with the 3-dimensional block technique. *Implant Dent* (2012) 21:444-8.
16. Araujo PP, Oliveira KP, Montenegro SC, Carreiro AF, Silva JS, Germano AR. Block allograft for reconstruction of alveolar bone ridge in implantology: a systematic review. *Implant Dent* (2013) 22:304-8.
17. Schlee M, Dehner J-F, Baukloh K, Happe A, Seitz O, Sader R. Esthetic outcome of implant-based reconstructions in augmented bone: comparison of autologous and allogeneic bone block grafting with the pink esthetic score (PES). *Head & Face Medicine* (2014) 10:21.
18. Monje A, Pikos MA, Chan HL, Suarez F, Gargallo-Albiol J, Hernandez-Alfaro F, Galindo-Moreno P, Wang HL. On the feasibility of utilizing allogeneic bone blocks for atrophic maxillary augmentation. *Biomed Res Int* (2014). Article ID 814578 doi:10.1155/2014/814578.
19. Leong DJ, Oh TJ, Benavides E, Al-Hezaimi K, Misch CE, Wang HL. Comparison between sandwich bone augmentation and allogenic block graft for vertical ridge augmentation in the posterior mandible. *Implant Dent* (2015) 24:4-12.
20. Laino L, Iezzi G, Piattelli A, Lo Muzio L, Cicciu M. Vertical ridge augmentation of the atrophic posterior mandible with sandwich technique: bone block from the chin area versus corticocancellous bone block allograft-clinical and histological prospective randomized controlled study. *Biomed Res Int* (2014). Article ID 982104 doi:10.1155/2014/982104.
21. Checchi V, Felice P. The inlay technique with an allograft block used for the rehabilitation of an atrophic posterior mandible: a case report. *Clinical Advances in Periodontics* (2015) 5:68-74.

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1. Eufinger H, Wehmoller M, Machtens E, Heuser L, Harders A, Kruse D. Reconstruction of craniofacial bone defects with individual alloplastic implants based on CAD/CAM manipulated CT-Data. *J. Craniomaxillofac. Surg.* (1995) 23:175-181.
2. Peters F, Groisman D, Davids R, Hanel T, Durr H, Klein M. Comparative study of patient individual implants from beta-tricalcium phosphate made by different techniques based on CT data. *Materialwiss. Werkstofftech.* (2006) 37:457-461.
3. Blume O, Back M, Born T, Smeets R, Jung O, Barbeck M. Treatment of a bilaterally severely resorbed posterior mandible due to early tooth loss by Guided Bone Regeneration using customized allogeneic bone blocks: A case report with 24 months follow-up data. *J Esthet Restor Dent* (2018), doi:10.1111/jerd.12388-
4. Klein M, Glatzer C. Individual CAD/CAM fabricated glass-bioceramic implants in reconstructive surgery of the bony orbital floor. *Plast Reconstr Surg* (2006) 117:565-570.
5. Hallermann W, Olsen S, Bardyn T, Taghizadeh F, Banic A, Iizuka T. A new method for computer-aided operation planning for extensive mandibular reconstruction. *Plast Reconstr Surg* (2006) 117:2431-7.
6. Khalyfa A, Vogt S, Weisser J, Grimm G, Rechtenbach A, Meyer W, Schnabelrauch M. Development of a new calcium phosphate powder-binder system for the 3D printing of patient specific implants. *J Mater Sci Mater Med* (2007) 18:909-16.
7. Binder WJ. Custom-designed facial implants. *Facial plastic surgery clinics of North America* (2008) 16:133-46, vii.
8. Cohen A, Laviv A, Berman P, Nashef R, Abu-Tair J. Mandibular reconstruction using stereolithographic 3-dimensional printing modeling technology. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* (2009) 108:661-6.
9. Guo LJ, Tian WD, Feng F, Long J, Li PF, Tang W. Reconstruction of orbital floor fractures comparison of individual prefabricated titanium implants and calvarial bone grafts. *Ann Plas Surg* (2009) 63:624-631.
10. Zhou L-B, Shang H-T, He L-S, Bo B, Liu G-C, Liu Y-P, Zhao J-L. Accurate reconstruction of discontinuous mandible using a reverse engineering/computer-aided design/rapid prototyping technique: a preliminary clinical study. *J Oral Maxillofac Surg* (2010) 68:2115-2121.
11. Lethaus B, Kessler P, Boeckman R, Poort LJ, Tolba R. Reconstruction of a maxillary defect with a fibula graft and titanium mesh using CAD/CAM techniques. *Head Face Med* (2010) 6:16.
12. Saldarriaga JFI, Vélez SC, Posada MDaC, Henao IEBB, Valencia MECaT. Design and manufacturing of a custom skull implant. *American Journal of Engineering & Applied Sciences* (2011) 4:169-174.
13. Brune JC, Hesselbarth U, Seifert P, Nowack D, Von Versen R, Smith MD, Seifert D. CT lesion model-based structural allografts: custom fabrication and clinical experience. *Transfusion Medicine and Hemotherapy* (2012) 39:395-404.
14. Ciocca L, Donati D, Fantini M, Landi E, Piattelli A, Iezzi G, Tampieri A, Spadari A, Romagnoli N, Scotti R. CAD-CAM-generated hydroxyapatite scaffold to replace the mandibular condyle in sheep: Preliminary results. *J Biomater Appl* (2013) 28:207-218.

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
Kroiss S. et al. Is socket healing conditioned by buccal plate thickness? A clinical and histologic study 4 months after mineralized human bone allografting. Clin Oral Implants Res (2014) 25:e120-6.

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Beck T.M. et al. Histologic analysis of healing after tooth extraction with ridge preservation using mineralized human bone allograft. J Periodontol (2010) 81:1765-1772.

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

 Bibliografie

1. Block MS, Finger I, Lytle R. Human mineralized bone in extraction sites before implant placement: Preliminary results. *J Am Dentn Assoc* (2002) 133:1631-1638.
2. Minichetti JC, D'amore JC, Hong AY, Cleveland DB. Human histologic analysis of mineralized bone allograft (Puros) placement before implant surgery. *J Oral Implantol* (2004) 30:74-82.
3. Wang HL, Kiyonobu K, Neiva RF. Socket augmentation: rationale and technique. *Implant Dent* (2004) 13:286-96.
4. Minichetti JC, D'amore JC, Hong AY. Three-year analysis of Tapered Screw-Vent implants placed into extraction sockets grafted with mineralized bone allograft. *J Oral Implantol* (2005) 31:283-93.
5. Jackson BJ, Morcos I. Socket grafting: a predictable technique for site preservation. *J Oral Implantol* (2007) 33:353-64.
6. Wang HL, Tsao YP. Mineralized bone allograft-plug socket augmentation: Rationale and technique. *Implant Dent* (2007) 16:33-41.
7. Danesh-Meyer M. Management of the extraction socket: site preservation prior to implant placement. *Australasian Dental Practice* (2008) 150-158.
8. Wang HL, Tsao YP. Histologic evaluation of socket augmentation with mineralized human allograft. *Int J Periodont Rest Dent* (2008) 28:231-7.
9. Fotek PD, Neiva RF, Wang HL. Comparison of dermal matrix and polytetrafluoroethylene membrane for socket bone augmentation: a clinical and histologic study. *J Periodontol* (2009) 80:776-85.
10. Beck TM, Mealey BL. Histologic analysis of healing after tooth extraction with ridge preservation using mineralized human bone allograft. *J Periodontol* (2010) 81:1765-72.
11. El Chaar ES. Soft tissue closure of grafted extraction sockets in the posterior maxilla: the rotated pedicle palatal connective tissue flap technique. *Implant Dent* (2010) 19:370-7.
12. Fickl S. Socket Preservation zum Zeitpunkt der Zahnextraktion. *Parodontologie Nachrichten* (2010) 1:14.
13. Kistler S. Praxisgerechte Techniken der Alveolar Ridge Preservation. *dentalspiegel* (2011) 5:22-26.
14. Tolstunov L, Chi J. Alveolar ridge augmentation: comparison of two socket graft materials in implant cases. *Compendium of Continuing Education in Dentistry* (2011) 32:E16-E124.
15. Fischer K, Jockel-Schneider Y, Bechtold M, Fickl S, Schlagenhauf U. „Socket preservation“ nach Zahnextraktion. *Der Freie Zahnarzt* (2011) 55:62-68.
16. Sterio TW, Katancik JA, Blanchard SB, Xenoudi P, Mealey BL. A prospective, multicenter study of bovine pericardium membrane with Allograft Cancellous Particles for localized alveolar ridge augmentation. *Int J Perio Rest Dent* (2013) 33:499-507.
17. Richter O. Tatort ästhetische Zone. *teamwork* (2013) 1:76-81.
18. Al-Hezaimi K, Rudek I, Al-Hamdan KS, Javed F, Nooh N, Wang HL. Efficacy of using a dual layer of membrane (dPTFE placed over collagen) for ridge preservation in fresh extraction sites: a micro-computed tomographic study in dogs. *Clin Oral Implants Res* (2013) 24:1152-7.
19. Spinato S, Galindo-Moreno P, Zaffe D, Bernardello F, Soardi CM. Is socket healing conditioned by buccal plate thickness? A clinical and histologic study 4 months after mineralized human bone allografting. *Clin Oral Implants Res* (2014) 25:e120-6.
20. Fu JH, Rios H, Al-Hezaimi K, Oh TJ, Benavides E, Wang HL. A randomized clinical trial evaluating the efficacy of the sandwich bone augmentation technique in increasing buccal bone thickness during implant placement. II. Tomographic, histologic, immunohistochemical, and RNA analyses. *Clin Oral Implants Res* (2014).
21. Fischer K, Fickl S. Alveolenmanagement: Eine Übersicht: Von Sofortimplantation bis Socket-Seal-Technik. *teamwork J CONT DENT EDUC* (2014) 6-14.

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Thorwarth M. et al. Evaluation of substitutes for bone: comparison of microradiographic and histological assessments. Br J Oral Max Surg (2007) 45:41-47.

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




Bibliografie

1. Block MS, Degen M. Horizontal ridge augmentation using human mineralized particulate bone: Preliminary results. *J Oral Maxillofac Surg* (2004) 62:67-72.
2. Minichetti JC, D'amore JC, Hong AY. Three-year analysis of Tapered Screw-Vent implants placed into extraction sockets grafted with mineralized bone allograft. *J Oral Implantol* (2005) 31:283-93.
3. Block MS. Horizontal ridge augmentation using particulate bone. *Atlas Oral Maxillofac Surg Clin North Am* (2006) 14:27-38.
4. Plöger M, Schau I. Rekonstruktion von Kieferkammdefekten mittels allogener Spongiosa. *Implantologie Zeitung* (2006) 10:16-18.
5. Park SH, Wang HL. Management of localized buccal dehiscence defect with allografts and acellular dermal matrix. *Int J Periodont Rest Dent* (2006) 26:589-95.
6. Schau I, Plöger M. Sofortimplantation und laterale Augmentation mit humanem Knochenersatzmaterial. *Dentale Implantologie & Parodontologie* (2006) 10:174-179.
7. Steigmann M. Barrieremembran zur Augmentation horizontaler Alveolarkammdefekte. *Zahnarzt & Praxis* (2007) 10:476-479.
8. Park SH, Lee KW, Oh TJ, Misch CE, Shotwell J, Wang HL. Effect of absorbable membranes on sandwich bone augmentation. *Clin Oral Implants Res* (2008) 19:32-41.
9. Le BT, Woo I. Alveolar cleft repair in adults using guided bone regeneration with mineralized allograft for dental implant site development: a report of 2 cases. *J Oral Maxillofac Surg* (2009) 67:1716-1722.
10. Lee A, Brown D, Wang HL. Sandwich bone augmentation for predictable horizontal bone augmentation. *Implant Dent* (2009) 18:282-90.
11. Burger BW. Use of ultrasound-activated resorbable poly-d-l-lactide pins (sonicpins) and foil panels (resorb-x) for horizontal bone augmentation of the maxillary and mandibular alveolar ridges. *J Oral Maxillofac Surg* (2010) 68:1656-1661.
12. Le B, Rohrer MD, Prasad HS. Screw "Tent-Pole" grafting technique for reconstruction of large vertical alveolar ridge defects using human mineralized allograft for implant site preparation. *J Oral Maxillofac Surg* (2010) 68:428-435.
13. Tolstunov L, Chi J. Alveolar ridge augmentation: comparison of two socket graft materials in implant cases. *Compendium Of Continuing Education In Dentistry* (2011) 32:E16-E124.
14. Fu JH, Wang HL. Horizontal bone augmentation: the decision tree. *Int J Periodontics Restorative Dent* (2011) 31:429-36.
15. Solakoglu Ö. Präimplantologische laterale Kieferkammaugmentation mit allogenen Knochenersatzmaterial: Eine Fallserie mit histologischer und histomorphometrischer Dokumentation. *Z Zahnärztl Impl* (2012) 27:24-32.
16. Steigmann M, Salama M, Wang H-L. Periosteal pocket flap for horizontal bone regeneration: a case series. *Int J Periodont Rest Dent* (2012) 32:311-320.
17. Fu J-H, Wang H-L. The sandwich bone augmentation technique. *Clinical Advances in Periodontics* (2012) 2:172-177.
18. Sterio TW, Katancik JA, Blanchard SB, Xenoudi P, Mealey BL. A prospective, multicenter study of bovine pericardium membrane with Allograft Cancellous Particles for localized alveolar ridge augmentation. *Int J Perio Rest Dent* (2013) 33:499-507.
19. Richter O. Tatort ästhetische Zone. *teamwork* (2013) 1:76-81.
20. Ronda M, Rebaudi A, Torelli L, Stacchi C. Expanded vs. dense polytetrafluoroethylene membranes in vertical ridge augmentation around dental implants: a prospective randomized controlled clinical trial. *Clin Oral Implants Res* (2014) 25:859-66.
21. Fu JH, Oh TJ, Benavides E, Rudek I, Wang HL. A randomized clinical trial evaluating the efficacy of the sandwich bone augmentation technique in increasing buccal bone thickness during implant placement surgery: I. Clinical and radiographic parameters. *Clin Oral Implants Res* (2014) 25:458-67.
22. Spinato S, Galindo-Moreno P, Zaffe D, Bernardello F, Soardi CM. Is socket healing conditioned by buccal plate thickness? A clinical and histologic study 4 months after mineralized human bone allografting. *Clin Oral Implants Res* (2014) 25:e120-6.

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Schmitt C.M. et al. Histological results after maxillary sinus augmentation with Straumann® BoneCeramic, Bio-Oss, Puros, and autologous bone. A randomized controlled clinical trial. Clin Oral Implants Res (2013) 24:576-585.

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
Froum S.J. et al. Comparison of mineralized cancellous bone allograft (Puros) and anorganic bovine bone matrix (Bio-Oss) for sinus augmentation: histomorphometry at 26 to 32 weeks after grafting. *Int J Periodont Rest* (2006) 26:543-551.

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
 Bibliografie

1. Tadic D, Beckmann F, Donath T, Epple M. Comparison of different methods for the preparation of porous bone substitution materials and structural investigations by synchrotron μ -computer tomography. *Materialwiss. Werkstofftech.* (2004) 35:240-244.
2. Thorwarth M, Schlegel KA, Srour S, Wiltfang J. Knöcherner Reparatur ossärer Defekte - Einsatz partikulärer homologer und boviner Knochenersatzmaterialien - erste tierexperimentelle Ergebnisse. *ZWR* (2004) 113:427-436.
3. Tadic D, Epple M. A thorough physicochemical characterisation of 14 calcium phosphate-based bone substitution materials in comparison to natural bone. *Biomaterials* (2004) 25:987-94.
4. Froum SJ, Tarnow DP, Wallace SS, Jalbout Z, Cho SC, Rohrer MD, Prasad HS. The use of a mineralized allograft for sinus augmentation: an interim histological case report from a prospective clinical study. *Compend Contin Educ Dent* (2005) 26:259-265.
5. Noubissi SS, Lozada JL, Boyne PJ, Rohrer MD, Clem D, Kim JS, Prasad H. Clinical, histologic, and histomorphometric evaluation of mineralized solvent-dehydrated bone allograft (Puros) in human maxillary sinus grafts. *J Oral Implantol* (2005) 31:171-9.
6. Schau I, Plöger M. Externer Sinuslift step-by-step: die laterale Window-Technik. *Dentale Implantologie & Parodontologie* (2005) 9:548-555.
7. Shin HI, Sohn DS. A method of sealing perforated sinus membrane and histologic finding of bone substitutes: a case report. *Implant Dent.* (2005) 14:328-33.
8. Froum SJ, Wallace SS, Elian N, Cho SC, Tarnow DP. Comparison of mineralized cancellous bone allograft (Puros) and anorganic bovine bone matrix (Bio-Oss) for sinus augmentation: histomorphometry at 26 to 32 weeks after grafting. *Int J Periodont Rest Dent* (2006) 26:543-51.
9. Gapski R, Neiva R, Oh TJ, Wang HL. Histologic analyses of human mineralized bone grafting material in sinus elevation procedures: a case series. *Int. J Periodont. Rest. Dent.* (2006) 26:59-69.
10. Whitesides LM, Radwan A, Sharawy M. Sinus floor augmentation using a composite graft of bone morphogenetic protein-2 and allogenic cancellous bone (Puros): case report. *J Oral Implantol.* (2006) 32:259-264.
11. Kneser U, Stangenberg L, Ohnolz J, Buettner O, Stern-Straeter J, Möbest D, Horch RE, Stark GB, Schaefer DJ. Evaluation of processed bovine cancellous bone matrix seeded with syngenic osteoblasts in a critical size calvarial defect rat model. *J Cell MolnMed* (2006) 10:695-707.
12. Aghaloo TL, Moy PK. Which hard tissue augmentation techniques are the most successful in furnishing bony support for implant placement? *Int J Oral Maxillofac Implants* (2007) 22 Suppl:49-70.
13. Thorwarth M, Wehrhan F, Srour S, Schultze-Mosgau S, Felszeghy E, Bader RD, Schlegel KA. Evaluation of substitutes for bone: comparison of microradiographic and histological assessments. *Br J Oral Maxillofac Surg* (2007) 45:41-7.
14. Baumgarten S. Allograft-Sinuslift mit XIVE-S plus Implantat. *Implant/Direkt* (2008) 3:6-9.
15. Minichetti JC, D'amore JC, Hong AY. Three-year analysis of tapered screw vent implants placed into maxillary sinuses grafted with mineralized bone allograft. *J Oral Implantol* (2008) 34:135-41.
16. Pikos MA. Maxillary sinus membrane repair: update on technique for large and complete perforations. *Implant Dent.* (2008) 17:24-31.
17. Meyer S, Floerkemeier T, Windhagen H. Histological osseointegration of Tutobone (R): first results in human. *Arch Orthop Traum Su* (2008) 128:539-544.
18. Tudor C, Srour S, Thorwarth M, Stockmann P, Neukam FW, Nkenke E, Schlegel KA, Felszeghy E. Bone regeneration in osseous defects - application of particulated human and bovine materials. *Oral Surg Oral Med O* (2008) 105:430-436.
19. Sohn DS, Lee JK, An KM, Shin HI. Histomorphometric Evaluation of Mineralized Cancellous Allograft in the Maxillary Sinus Augmentation: A 4 Case Report *Implant Dent* (2009) 18:172-181.

20. Baumgarten S. Sinus grafting with mineralized allograft and staged implant placement. *Quintessence Int.* (2010) 41:197-202.
21. Noubissi SS, Kleinman AJ. A comparative clinical, histological and histomorphometric study of mineralized allograft and xenograft materials in the treatment of atrophic maxillary sinuses. *Journal of Implant and Advanced Clinical Dentistry* (2010) 2:47-60.
22. Zimmer Dental. Der Sinuslift: Risiken und Möglichkeiten eines etablierten Routineeingriffs. *ZWR* (2010) 119:450,450.
23. Annibali S, Cristalli MP, La Monaca G, Bignozzi I, Scarano A, Corrado R, Lo Muzio L. Human maxillary sinuses augmented with mineralized, solvent-dehydrated bone allograft: a longitudinal case series. *Implant Dent* (2011) 20:445-54.
24. Soardi CM, Spinato S, Zaffe D, Wang HL. Atrophic maxillary floor augmentation by mineralized human bone allograft in sinuses of different size: an histologic and histomorphometric analysis. *Clin Oral Implants Res* (2011) 22:560-6.
25. Bernardello F, Righi D, Cosci F, Bozzoli P, Soardi Carlo M, Spinato S. Crestal sinus lift with sequential drills and simultaneous implant placement in sites with <5 mm of native bone: a multicenter retrospective study. *Implant Dent* (2011) 20:439-444 10.1097/ID.0b013e3182342052.
26. Soardi CM, Wang H-L. New crestal approach for lifting sinus in the extremely atrophic upper maxillae. *Clinical Advances in Periodontics* (2012) 2:179-185.
27. Ayobian-Markazi N, Fourootan T, Kharazifar MJ. Comparison of cell viability and morphology of a human osteoblast-like cell line (SaOS-2) seeded on various bone substitute materials: An in vitro study. *Dental research journal* (2012) 9:86-92.
28. Schmitt CM, Doering H, Schmidt T, Lutz R, Neukam FW, Schlegel KA. Histological results after maxillary sinus augmentation with Straumann(R) BoneCeramic, Bio-Oss(R), Puros(R), and autologous bone. A randomized controlled clinical trial. *Clin Oral Implants Res* (2013) 24:576-85.
29. Berberi A, Nader N. Subantral augmentation with mineralized cortical bone allograft material: clinical, histological, and histomorphometric analyses and graft volume assessments. *Implant Dent* (2016) 25:353-60.
30. Soardi CM, Clozza E, Turco G, Biasotto M, Engebretson SP, Wang H-L, Zaffe D. Microradiography and microcomputed tomography comparative analysis in human bone cores harvested after maxillary sinus augmentation: a pilot study. *Clin Oral Implants Res* (2014) 25:1161-1168.

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Browning E.S. et al. Evaluation of a mineralized cancellous bone allograft for the treatment of periodontal osseous defects: 6-month surgical reentry. Int J Periodont Rest Dent (2009) 29:41-47.

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
Bibliografie

1. Yukna RA, Vastardis S, Mayer ET. Puros vs. Demineralized freeze-dried bone in human periodontal osseous defects. in IADR/AADR/CADR 83rd General Session. 2005. Baltimore.
2. Tsao YP, Neiva R, Al-Shammari K, Oh TJ, Wang HL. Effects of a mineralized human cancellous bone allograft in regeneration of mandibular Class II furcation defects. *J Periodontol* (2006) 77:416-25.
3. Vastardis S, Yukna RA. Evaluation of allogeneic bone graft substitute for treatment of periodontal osseous defects: 6-month clinical results. *Compend Contin Educ Dent* (2006) 27:38-44.
4. Wang HL, Boyapati L. "PASS" principles for predictable bone regeneration. *Implant Dent* (2006) 15:8-17.
5. Browning ES, Mealey BL, Mellonig JT. Evaluation of a mineralized cancellous bone allograft for the treatment of periodontal osseous defects: 6-month surgical reentry. *Int J Periodont Rest Dent* (2009) 29:41-47.
6. Koylass JM, Valderrama P, Mellonig JT. Histologic evaluation of an allogeneic mineralized bone matrix in the treatment of periodontal osseous defects. *Int J Periodontics Restorative Dent* (2012) 32:405-11.
7. Engler-Hamm D, Fickl S. Langzeiterfolg der parodontalen Regeneration. *Quintessenz* (2015) 66:761-767.

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1. Daten in Akten bei RTI Biologics, Inc.
ZimVie Germany GmbH
ZimVie Austria GmbH
Biomet 3i Schweiz GmbH