




Product Catalog
2023



GenTek™

Genuine Restorative Components

Certain®, External Hex, and Low Profile





GenTek Open Connectivity

For even more flexibility.

Genuine

GenTek Pre-Abutment Blanks are compatible with Medentika® holders (CAM update required).



NEW

GenTek TiBases


Achieve efficiency to your CAD/CAM restoration workflow with a singular TiBase that provides flexibility of channel angulation and adjustable post height options.

exocad

RealGUIDE™
UNIVERSAL OPEN SYSTEM

dental wings

3shape



Free

Libraries for the GenTek restorative components are available for Zfx CAD, RealGUIDE™, 3Shape, exocad, and Dental Wings CAD software.

Free download at zfx-dental.com



SCAN & INFORM
More at
www.zfx-dental.com








Table of Contents

NEW: GenTek TiBases	4
GenTek Pre-Abutment Blanks	6
GenTek TiBases – Milling Block Compatible	8
GenTek TiBases – ZimVie Multi-Unit Abutment Compatible	12
GenTek Scanbodies	14
GenTek Digital Scan Analogs	17

Legend

∅ Restorative Platform (mm)
Color code for Certain
and External Hex

	3.4 mm
	4.0–6.0 mm
	5.0 mm
	6.0 mm

Placement Tool for GenTek Digital Scan Analogs Compatibility Guide

Placement Tool	Laboratory Analog
M 14	Low Profile Abutment Lab Analogs
M 16	Certain Lab Analogs
M 20	External Hex Lab Analogs



GenTek Calibration Tool for Milling Machines

Description	Part No.
GenTek Calibration Abutment Blank	ZFX18-CALIB



GenTek Pre-Abutment Blank Holder

Description	Part No.
Zfx Pre-Abutment Blank Holder, Standard	ZFX02PMSH
Zfx Pre-Abutment Blank Holder, Multi-Blank for Zfx Inhouse5x	ZFX02002069
Zfx Pre-Abutment Blank Adaptor for Glass Ceramic Holder	ZFX02PMA




Milling Block Compatibility Guide¹

ZimVie Restorative Platform Connection	Block Interface Size
Certain	S / L ²
External Hex	L



¹ Not applicable for the GenTek Angulated Screw Channel TiBase

² Restorative Platform 3.4 mm is compatible with interface size S
Restorative Platform 4.1 and 5.0 mm are compatible with interface size L



NEW: GenTek TiBases

ZimVie Dental Certain Restorative Platform Connection

TiBase Engaging

incl. Gold-Tite Hexalobular Screw and Gluing Screw



Post Height ¹ (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.	
7	3.4	0.3	ZFX11ZB-CE34ASZ03E	
		1.3	ZFX11ZB-CE34AS13E	
		2.6	ZFX11ZB-CE34AS26E	
	4.1	4.1	0.3	ZFX11ZB-CE41ASZ03E
			1.3	ZFX11ZB-CE41AS13E
			2.6	ZFX11ZB-CE41AS26E
	5.0	5.0	0.3	ZFX11ZB-CE50ASZ03E
			1.3	ZFX11ZB-CE50AS13E
			2.6	ZFX11ZB-CE50AS26E

Gold-Tite® Hexalobular Screw

Hexalobular Try-In Screw



Part No.		Part No.	
ZFX09-ZB-CE-HLGTE	20 Ncm	ZFX09-ZB-CE-HLTISE	Hand-tighten

TiBase Non-Engaging

incl. Gold-Tite Hexalobular Screw and Gluing Screw



Post Height ¹ (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.
7	3.4	0.3	ZFX11ZB-CE34ASZ03N
		1.3	ZFX11ZB-CE34AS13N
	4.1	0.3	ZFX11ZB-CE41ASZ03N
		1.3	ZFX11ZB-CE41AS13N
	5.0	0.3	ZFX11ZB-CE50ASZ03N
		1.3	ZFX11ZB-CE50AS13N

Gold-Tite Hexalobular Screw

Hexalobular Try-In Screw



Part No.		Part No.	
ZFX09-ZB-CE-HLGTN	20 Ncm	ZFX09-ZB-CE-HLTISN	Hand-tighten

¹ Various groove heights at 4.5, 5.7 mm

Tools for New GenTek TiBases

Hexalobular Screwdriver for Hexalobular Retaining Screw



Length (mm)	Part No.
21	ZFX02HLD21
28	ZFX02HLD28
36	ZFX02HLD36

Low Torque Indicating Ratchet Wrench



Description	Part No.
Low Torque Indicating Ratchet Wrench	L-TIRW
L-TIRW Standard ISO 1797 Adapter	C9980

Certain QuickSeat® Activator Tool



Part No.
IQSA01

Low Torque Indicating Ratchet Wrench Kit¹ incl. Adapter, Tray, and Driver Tips



Description	Part No.
<ul style="list-style-type: none"> L-TIRW: Low Torque Indicating Ratchet Wrench C9980: L-TIRW Standard ISO 1797 Adapter^{2,3} RTI2035TR: Torque Indicating Ratchet Wrench Tray RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24mm RASQ3N: Narrow Right Angle Square Driver Tip, 24mm 	L-TIRWK

¹ Components also sold separately
² Short ISO 1797 adapter also available (C9981)
³ C9980 can be ordered separately

GenTek Pre-Abutment Blanks

ZimVie Dental Certain Restorative Platform Connection

Titanium Pre-milled Abutment Blank *incl. Gold-Tite Screw*



Restorative Platform (mm)	Part No.
3.4	ZFX18-ZB-CE-34
4.1, 5.0, 6.0	ZFX18-ZB-CE-4160

Gold-Tite Screw



Part No.	Recommended Torque
ZFX09-ZB-CE-HGTS	20 Ncm

Try-In Screw



Part No.	Recommended Torque
ZFX09-ZB-CE-HTIS	Hand-tighten

Hex Driver, Driver Tip *for Gold-Tite and Try-In Screws*



Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

Low Torque Indicating Ratchet Wrench Kit¹ *incl. Adapter, Tray, and Driver Tips*



Description	Part No.
<ul style="list-style-type: none"> L-TIRW: Low Torque Indicating Ratchet Wrench C9980: L-TIRW Standard ISO 1797 Adapter^{2/3} RTI2035TR: Torque Indicating Ratchet Wrench Tray RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24mm RASQ3N: Narrow Right Angle Square Driver Tip, 24mm 	L-TIRWK

Certain QuickSeat Activator Tool



Part No.
IQSA01

ZimVie Dental External Hex Restorative Platform Connection

Titanium Pre-milled Abutment Blank
incl. Gold-Tite Screw



Restorative Platform (mm)	Part No.
3.4	ZFX18-ZB-EX-34
4.1, 5.0, 6.0	ZFX18-ZB-EX-4160

Gold-Tite Screw

Try-In Screw



Recommended Torque



Recommended Torque

Gold-Tite Screw		Try-In Screw	
Part No.	Recommended Torque	Part No.	Recommended Torque
ZFX09-ZB-EX-GTS	32 Ncm	ZFX09-ZB-EX-TIS	Hand-tighten

Square Driver, Driver Tip
for Gold-Tite Screw



Description	Part No.
Narrow Posterior Square Driver, 17 mmL	PSQD0N
Narrow Standard Slotted Square Driver, 24 mmL	PSQD1N
Narrow Right Angle Square Driver Tip, 24 mmL	RASQ3N
Narrow Right Angle Square Driver Tip, 30 mmL	RASH8N

Hex Driver, Driver Tip
for Try-In Screw



Only for External Hex Try-In Screw

Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

Low Torque Indicating Ratchet Wrench Kit¹
incl. Adapter, Tray, and Driver Tips



Description	Part No.
<ul style="list-style-type: none"> L-TIRW: Low Torque Indicating Ratchet Wrench C9980: L-TIRW Standard ISO 1797 Adapter^{2,3} RTI2035TR: Torque Indicating Ratchet Wrench Tray RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24 mm RASQ3N: Narrow Right Angle Square Driver Tip, 24 mm 	L-TIRWK

¹Components also sold separately
²Short ISO 1797 adapter also available (C9981)
³C9980 can be ordered separately

GenTek TiBases Milling Block Compatible

ZimVie Dental Certain Restorative Platform Connection

Single-Unit Restorative Components

TiBase Engaging¹
incl. Gold-Tite and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-CE-3447-ES
	4.1	ZFX11-ZB-CE-4147-EL
	5.0	ZFX11-ZB-CE-5047-EL

Gold-Tite Screw

Try-In Screw



Part No.		Part No.	
ZFX09-ZB-CE-HGTS	20 Ncm	ZFX09-ZB-CE-HTIS	Hand-tighten

Multi-Unit Restorative Components

TiBase Non-Engaging¹
incl. Gold-Tite and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-CE-3447-NES
	4.1	ZFX11-ZB-CE-4147-NEL
	5.0	ZFX11-ZB-CE-5047-NEL

Gold-Tite Screw

Try-In Screw



Part No.		Part No.	
ZFX09-ZB-CE-LGTS	20 Ncm	ZFX09-ZB-CE-NHTIS	Hand-tighten

Hex Driver, Driver Tip
for Gold-Tite and Try-In Screws



Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

Low Torque Indicating Ratchet Wrench Kit²
incl. Adapter, Tray, and Driver Tips



Certain QuickSeat[®] Activator Tool



Description	Part No.	Part No.
<ul style="list-style-type: none"> L-TIRW: Low Torque Indicating Ratchet Wrench C9980: L-TIRW Standard ISO 1797 Adapter^{3,4} RTI2035TR: Torque Indicating Ratchet Wrench Tray RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24mm RASQ3N: Narrow Right Angle Square Driver Tip, 24mm 	L-TIRWK	IQSA01


¹4 mm groove indicating minimal allowed post height

²Not applicable for the GenTek Angulated Screw Channel TiBase

³Components also sold separately

⁴Short ISO 1797 adapter also available (C9981)

⁵C9980 can be ordered separately



GenTek TiBases Milling Block Compatible

ZimVie Dental External Hex Restorative Platform Connection

Single-Unit Restorative Components

TiBase Engaging¹
incl. Gold-Tite and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-EX-3447-EL
	4.1	ZFX11-ZB-EX-4147-EL
	5.0	ZFX11-ZB-EX-5047-EL

Gold-Tite Screw



Recommended Torque

Part No.	
ZFX09-ZB-EX-GTS	32 Ncm

Try-In Screw



Recommended Torque

Part No.	
ZFX09-ZB-EX-TIS	Hand-tighten

Multi-Unit Restorative Components

TiBase Non-Engaging¹
incl. Gold-Tite and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-EX-3447-NEL
	4.1	ZFX11-ZB-EX-4147-NEL
	5.0	ZFX11-ZB-EX-5047-NEL

Gold-Tite Screw



Recommended Torque

Part No.	
ZFX09-ZB-EX-GTS	32 Ncm

Try-In Screw



Recommended Torque

Part No.	
ZFX09-ZB-EX-TIS	Hand-tighten

Square Driver, Driver Tip
for Gold-Tite Screw



Hex Driver, Driver Tip
for Try-In Screw



Only for
External Hex
Try-In Screw

Description	Part No.
Narrow Posterior Square Driver, 17 mmL	PSQD0N
Narrow Standard Slotted Square Driver, 24 mmL	PSQD1N
Narrow Right Angle Square Driver Tip, 24 mmL	RASQ3N
Narrow Right Angle Square Driver Tip, 30 mmL	RASH8N

Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

Low Torque Indicating Ratchet Wrench Kit²
incl. Adapter, Tray, and Driver Tips



Description	Part No.
<ul style="list-style-type: none"> L-TIRW: Low Torque Indicating Ratchet Wrench C9980: L-TIRW Standard ISO 1797 Adapter^{3,4} RTI2035TR: Torque Indicating Ratchet Wrench Tray RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24 mm RASQ3N: Narrow Right Angle Square Driver Tip, 24 mm 	L-TIRWK

¹4 mm groove indicating minimal allowed post height

²Not applicable for the GenTek Angulated Screw Channel TiBase

³Components also sold separately

⁴Short ISO 1797 adapter also available (C9981)

⁵C9980 can be ordered separately

GenTek TiBases ZimVie Multi-Unit Abutment Compatible

ZimVie Dental Low Profile Abutment

Multi-Unit Restorative Components

TiBase Non-Engaging¹
incl. Gold-Tite and Gluing Screws




Height (mm)	Restorative Platform (mm)	Part No.
5.0	4.8	ZFX11-ZB-LP-48

Gold-Tite Screw

Try-in Screw



Part No.		Part No.
ZFX09-ZB-LP-GTS	10Ncm	ZFX09-ZB-LP-RS
		Hand-tighten



Hex Driver, Driver Tip
for Gold-Tite Screws and Try-in Screws



Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

Low Torque Indicating Ratchet Wrench Kit²
incl. Adapter, Tray, and Driver Tips



Description	Part No.
<ul style="list-style-type: none"> L-TIRW: Low Torque Indicating Ratchet Wrench C9980: L-TIRW Standard ISO 1797 Adapter^{3,4} RTI2035TR: Torque Indicating Ratchet Wrench Tray RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24 mm RASQ3N: Narrow Right Angle Square Driver Tip, 24 mm 	L-TIRWK

¹Compatible with non-hexed low profile abutment only
Not applicable for the GenTek Angulated Screw Channel TiBase
²Components also sold separately
³Short ISO 1797 adapter also available (C9981)
⁴C9980 can be ordered separately



ZimVie Dental Certain Restorative Platform Connection

Desktop Scanbody incl. Screw



Restorative Platform (mm)	Part No.	
3.4	ZFX05D-ZB-CE-34	Hand-tighten ¹
4.1, 5.0, 6.0	ZFX05D-ZB-CE-4160	
5.0	ZFX05D-ZB-CE-50 ²	
6.0	ZFX05D-ZB-CE-60 ²	

Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

Intraoral Scanbodies incl. Screws



Height (mm)	Restorative Platform (mm)	Part No.	
4.0 & 7.0	3.4	ZFX05I-ZB-CE-34	Hand-tighten ¹
	4.1, 5.0, 6.0	ZFX05I-ZB-CE-4160	

Hex Driver for Intraoral Scanbodies



Description	Part No.
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25

Certain QuickSeat® Activator Tool



Part No.
IQSA01

ZimVie Dental External Hex Restorative Platform Connection

Desktop Scanbody incl. Screw



Recommended Torque


Restorative Platform (mm)	Part No.	Hand-tighten ¹
3.4	ZFX05D-ZB-EX-34	
4.1, 5.0, 6.0	ZFX05D-ZB-EX-4160	
5.0	ZFX05D-ZB-EX-50 ²	
6.0	ZFX05D-ZB-EX-60 ²	

Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

Intraoral Scanbodies incl. Screws



Recommended Torque

Height (mm)	Restorative Platform (mm)	Part No.	Hand-tighten ¹
4.0 & 7.0	3.4	ZFX05I-ZB-EX-34	
	4.1, 5.0, 6.0	ZFX05I-ZB-EX-4160	

Hex Driver for Intraoral Scanbodies




Description	Part No.
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25

¹Minimum of 5 Ncm; Maximum of 15 Ncm.
Refer to Instruction Manual for more information.
²Only for use with Zfx Evolution scanners



ZimVie Dental Low Profile Abutment

Desktop Scanbody *incl. Screw*



↑
Recommended
Torque

Restorative Platform (mm)	Part No.		
4.8	ZFX05D-ZB-LP-48		Hand-tighten ¹

Torque Wrench *for Desktop Scanbody*



Part No.
ZFX02002067

Intraoral Scanbodies Engaging *incl. Screws*



↑
Recommended
Torque

Height (mm)	Restorative Platform (mm)	Part No.		
4.0	4.8	ZFX05I-ZB-LP-48		Hand-tighten ¹

Hex Driver, Latch-Lock Hex Driver *for Intraoral Scanbodies*



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short, 1.25 mmD, 23 mmL <i>(may require reduction hand piece)</i>	HX1.25D
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long, 1.25 mmD, 26 mmL <i>(may require reduction hand piece)</i>	HXL1.25D

Intraoral Scanbodies Non-Engaging *incl. Screws*



↑
Recommended
Torque

Height (mm)	Restorative Platform (mm)	Part No.		
4.0	4.8	ZFX05I-ZB-LP-48-NE		Hand-tighten ¹

¹ Minimum of 5 Ncm; Maximum of 15 Ncm.
Refer to Instruction Manual for more information.

GenTek Digital Scan Analogs

ZimVie Dental Certain Restorative Platform Connection

Digital Scan Analog¹



Restorative Platform (mm)	Part No.
3.4	ZFX15-ZB-CE-34
4.1	ZFX15-ZB-CE-41
5.0	ZFX15-ZB-CE-50
6.0	ZFX15-ZB-CE-60

Placement Tool for Digital Scan Analogs



Part No.
ZFX02-M16

ZimVie Dental External Hex Restorative Platform Connection

Digital Scan Analog¹



Restorative Platform (mm)	Part No.
3.4	ZFX15-ZB-EX-34
4.1	ZFX15-ZB-EX-41
5.0	ZFX15-ZB-EX-50
6.0	ZFX15-ZB-EX-60

Placement Tool for Digital Scan Analogs



Part No.
ZFX02-M20

ZimVie Dental Low Profile Abutment

Digital Scan Analog¹



Restorative Platform (mm)	Part No.
4.8	ZFX15-ZB-LP-48

Placement Tool for Digital Scan Analogs



Part No.
ZFX02-M14

¹Compatible with existing Certain / External Hex / Eztec / TSV/TM / Low Profile / Tapered Abutment Impression Copings

To order genuine ZimVie Dental Components
call your local representative or send an email to
dentalCS@ZimVie.com

ZimVie Dental Global Headquarters

4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272
www.ZimVie.com/dental



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. TiBases (Class IIb), Scanbodies (Class I), Blanks (Class IIb), Tools (Class I), Screws (Class IIb) and Analogs (non medical devices). ZV0217 REV B 05/23 ©2023 ZimVie. All rights reserved.

