

# RealGUIDE™ Z3D Guided Surgery

Designed for ZimVie's TSV® Trabecular Metal™ and Eztetic® Dental Implants





## RealGUIDE™ Z3D GUIDED SURGERY KIT

The new RealGUIDE Z3D Guided Surgery Kit contains instruments needed for the planned surgery. It is designed to work with ZimVie's **TSV, Trabecular Metal, and Eztetic** Implant Systems.

The kit information is available in the RealGUIDE Software Suite, an open system for diagnosis and advanced surgical guide design that will support you with customized drilling reports to be well-prepared for dental implant therapy.



RealGUIDE Software Suite  
with PLAN, APP, GUIDE Modules

- Fully guided placement of ZimVie Dental Eztetic 3.1mm Implant and TSV/TM Implants in sizes 3.7mm, 4.1mm, and 4.7mm.
- Color-coded for ease-of-use.
- Compact kit, which is useful for space-constrained dental offices.
- No handles or keys are necessary while preparing the osteotomy.
- DLC (Diamond-like Carbon) coating on the drills provides wear resistance characteristics .
- Includes implant mounts hex feature to align with the hex feature on the corresponding guide sleeves.
- Staging pedestal for holding the implant vial to remove the Fixture Mount Transfer (FMT) of the TSV implants is available. This enables the TSV implants to be placed fully guided through the surgical guides.
- Both soft and dense bone protocols are available.

Product	Part No.
RealGUIDE Z3D Surgical Kit	3DM0070Z



MORE INFORMATION  
ABOUT THE RealGUIDE Z3D  
GUIDED SURGERY KIT AT  
[ZIMVIE.COM](http://ZIMVIE.COM)

For more information, visit [ZimVie.com](http://ZimVie.com)

ZimVie  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
1-800-342-5454

 **ZimVie**  
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Holdings, Inc., and distributed and marketed by ZimVie Dental (and, in the case of distribution and marketing, its authorized marketing partners). The RealGUIDE Z3D is manufactured by 3Diemme, a subsidiary of ZimVie. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1217 REV A 05/23 ©2023 ZimVie. All rights reserved.

