


Stainless  
steel construct

# Dental Surgical Wash Tray



## Compartmentalized

Houses small, reusable instruments in two removable compartments, one with grommets to secure the latch-lock end of drills and tools, and a second utility compartment to hold other small instruments.

## Compatibility

The porous, stainless steel construct allows for detergent and water passage during the cleaning and rinsing processes in automated cleaners and thermo-disinfectors.

## Universal application

Holds small instruments for a variety of implant systems provided by ZimVie Dental.

 **ZimVie**

# Dental Surgical Wash Tray

The Dental Surgical Wash Tray is a stainless steel construct which holds small, reusable dental implantology instruments during the automated wash and disinfection process when carried out by a thermo-disinfector.


## Compartmentalized

- Wash tray contains two removable components
- Grommet compartment may hold up to 50 instruments with ISO-standard latch-lock shank
- Utility compartment with lid holds small, non-latch-lock instruments



## Compatibility

- Stainless steel construct compatible with enzymatic and alkaline solutions used by automated cleaners/thermo-disinfectors
- Porous design allows detergent and water passage during cleaning and rinsing process



## Universal

- Houses small implantology instruments for a variety of ZimVie dental implant systems during the automated cleaning process



## Ordering Information

Item Number	Description
ZBDWT01	Dental Surgical Wash Tray

For more information, visit [ZimVie.com](https://www.zimvie.com)

**ZimVie**  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
1-800-342-5454  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272  
[dentalCS@ZimVie.com](mailto:dentalCS@ZimVie.com)

 **ZimVie**  
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV2036 REV A 12/23 ©2023 ZimVie. All rights reserved.

