

GemLock® Delivery System



Performance that delivers

Secure Fit

Touchless delivery tools offer a reliable engagement with implants, healing abutments, and screws

Visual Indicators

Surgical drivers for TSX®, Trabecular Metal®, and Tapered Screw-Vent® Implants now include hex- and depth-indicating geometry and a longer length option

Expanded Offering

GemLock retention is available on restorative drivers for Encode® Healing Abutments and retaining screws for the Certain® connection

Consistent Performance

Durable synthetic rubies enhance driver tip longevity and ensure delivery tools remain consistent over multiple uses

Ordering Information

Placement Driver Tips for Certain Connection Implants

Description	Mode of Use	Item No.
Internal Connection Universal Placement Driver Tip, Short	Handpiece	IIPDTUS
Internal Connection Universal Placement Driver Tip, Long	Handpiece	IIPDTUL



Restorative Driver Tips, 1.2 mmD (0.048"), for Certain Healing Abutments and Abutment Screws

Description	Mode of Use	Item No.
Latch-lock Hex Driver with GemLock Retention, 1.2 mm (D), 24 mm (L)	Handpiece	IHX1.2D
Latch-lock Hex Driver with GemLock Retention, 1.2 mm (D), 30 mm (L)	Handpiece	IXL1.2D



Hex Drivers, 2.1 mmD, for 3.1 mmD TSX and Eztetic® Implants

Description	Mode of Use	Item No.
Latch-lock Hex Driver with GemLock Retention, 2.1 mm (D), Short	Handpiece	CHD2.1
Latch-lock Hex Driver with GemLock Retention, 2.1 mm (D), Long	Handpiece	CHDL2.1
Ratchet Hex Driver with GemLock Retention, 2.1 mm (D), Short	ZTIRW, RSR, SSHS	CHR2.1
Ratchet Hex Driver with GemLock Retention, 2.1 mm (D), Long	ZTIRW, RSR, SSHS	CHRL2.1



Hex Drivers, 2.5 mmD, for TSX, Trabecular Metal, and Tapered Screw-Vent Implants

Description	Mode of Use	Item No.
Latch-lock Hex Driver with GemLock Retention, 2.5 mm (D), 21 mm (L)	Handpiece	GHX2.5D
Latch-lock Hex Driver with GemLock Retention, 2.5 mm (D), 30 mm (L)	Handpiece	GHX2.5DL
Ratchet Hex Driver with GemLock Retention, 2.5 mm (D), 17 mm (L)	ZTIRW, RSR, SSHS	GHX2.5
Ratchet Hex Driver with GemLock Retention, 2.5 mm (D), 28 mm (L)	ZTIRW, RSR, SSHS	GHX2.5L



Restorative Drivers, 1.25 mmD, for TSX, Trabecular Metal, TSV, and Eztetic Implant Systems

Description	Mode of Use	Item No.
Hex Driver, Short, with GemLock Retention, 1.25 mm (D), 22 mm (L)	Hand	HXGR1.25
Hex Driver, Long, with GemLock Retention, 1.25 mm (D), 30 mm (L)	Hand	HXLGR1.25



For more information, visit ZimVie.com

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. ZTIRW High Torque Indicating Ratchet Wrench and L-TIRW Low Torque Indicating Ratchet Wrench are manufactured by Elos Medtech Pinol A/S. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV2105 REV A 02/24 ©2024 ZimVie. All rights reserved.

