

# Azure™ Tapered Bone Level Guided Surgical Kit



## Soaring to new heights in Guided Surgery

The Azure™ Tapered Bone Level (TBL) Guided Surgical Kit offers fully guided implant placement, designed to confidently craft clinical success placing all Azure TBL Implants.

## Productivity and Performance

### Immediacy and Stability

- Achieve primary stability for immediacy in all bone types for single to full-arch cases
- > 35 Ncm of Insertion Torque, even in very soft type IV bone with Ø3.3 mm implants<sup>1</sup>

### Clinical Confidence Creates Aesthetic Results

- Obtain maximum flexibility for beautiful patient outcomes
- Enabled by two guide sleeve diameters and two prolongations
- Now a guided solution for narrow anterior spaces and limited posterior cases
- Grow your practice with tailored patient-specific treatment plans
- DLC (Diamond-like Carbon) drills increase surface hardness and reduce frictional forces for precise osteotomies

### Cost-Effective, Easy-to-Use

- No need for cumbersome reduction handles
- Compact and easy-to-follow intuitive surgical kit
- Minimized number of drills and instruments, yet a comprehensive system to treat specific patient financial demographics

<sup>1</sup> Data on file with ZimVie. Data from testing on #20 type soft bone block, averages per diameter. Preclinical testing is not necessarily indicative of clinical results.

## Ordering Information

	Part Number	Description
<b>Guided</b>	AGTBL-GSURKIT*	Azure Tapered Bone Level Guided Surgical Kit, Complete
<b>Non-Guided</b>	ASTBL-SURKIT*	Azure Tapered Bone Level Implant Surgical Kit, Complete



\*Part number AGTBL-GSURKIT includes the tray, drills, and all instrumentation for guided surgery. Both the AGTBL-GSURKIT and ASTBL-SURKIT part numbers are for ordering purposes only. Individual kit parts are packaged separately.

Azure TBL Implants**	Ø3.3 mm (P)	Ø4.1 mm (P)	Ø4.8 mm (P)
	NC	RC	RC
<i>Implant-Lengths</i>			
8 mmL	AITBL-IMP3308	AITBL-IMP4108	AITBL-IMP4808
10 mmL	AITBL-IMP3310	AITBL-IMP4110	AITBL-IMP4810
12 mmL	AITBL-IMP3312	AITBL-IMP4112	AITBL-IMP4812
14 mmL	AITBL-IMP3314	AITBL-IMP4114	AITBL-IMP4814
16 mmL	AITBL-IMP3316	AITBL-IMP4116	AITBL-IMP4816



\*\* Cover Screw and Healing Abutments sold separately

Cover Screws	Ø3.3 mm (P)	Ø4.1–4.8mm (P)
	NC	RC
	ASTBL-CSNC33NS	ASTBL-CSRC4NS



Healing Abutments	Ø3.3 mm (P)	Ø4.1–4.8mm (P)
	NC	RC
<i>Cuff Height</i>		
3 mm	ARSBL-HA333NS	ARSBL-HA413NS
4 mm	ARSBL-HA334NS	ARSBL-HA414NS
5 mm	ARSBL-HA335NS	ARSBL-HA415NS



Please reference the Azure Tapered Bone Level Catalog (ZVINST0201) for complete ordering information.



Learn More

For more information, visit [azuredental.com](http://azuredental.com)

**azure** | A ZimVie Brand

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. and its affiliates; and all products are manufactured by one or more of the dental subsidiaries of ZimVie US Corp LLC (Biomet 3i, LLC, etc.) and distributed by ZimVie US Corp LLC and its authorized marketing partners. Azure Tapered Bone Level Implants, Restorative Components and ASTBL-SURKIT instruments are manufactured by Terrats Medical S.L. In the AGTBL-GSURKIT, instrumentation is manufactured by BIOTEC S.r.l., excluding the ASTBL-TIR, ASTBL-RA, ASTBL-IDA, ARA-DT20TC06, and ARA-DT30C06 manufactured by Terrats Medical S.L. EU Medical Device Classification: Azure TBL (Iib), Cover Screws (Iib), Restorative Components (IIa). The intended purpose of Azure TBL Implants is for the functional and esthetic oral rehabilitation of the upper or lower jaw. Healing Abutments are intended for the formation of a soft tissue sulcus. Cover Screws are to be temporarily connected to a dental implant to protect the implant connection during bone healing. The products shown in this material are not eligible for reimbursement in France. For other countries, please check the local laws and regulations. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV2868-GLBL REV B 10/25 ©2025 ZimVie.



The NB number applicable to each device can be found on the product label, if applicable.