



# Azure™ Tapered Bone Level Surgical Kit

## Usage Chart



Use a dry erase marker to track the total number of osteotomy preparations that are performed with each drill. See RILM (Reusable Instrument Lifespan Manual) for guidance on when to replace the drills.

\* Single use drills should be disposed of after each procedure.



Part #	1203G	ASTBL-ND1633	ASTBL-ND1641	ASTBL-PID2229	ASTBL-PID2233	ASTBL-PID2241	ASTBL-D2829	ASTBL-D2833	ASTBL-D2841	ASTBL-D3529	ASTBL-D3533	ASTBL-D3541	ASTBL-D4229	ASTBL-D4233	ASTBL-D4241	ASTBL-PRD3325	ASTBL-PRD3333	ASTBL-PRD4125	ASTBL-PRD4133	ASTBL-PRD4825	ASTBL-PRD4833	ASTBL-IDNC3328	ASTBL-IDNC3337	ASTBL-IDRC428	ASTBL-IDRC437
Product	Round Bur	Needle Drill 1.6 mmD x 33 mmL	Needle Drill 1.6 mmD x 41 mmL	Pilot Drill 2.2 mmD x 29 mmL	Pilot Drill 2.2 mmD x 33 mmL	Pilot Drill 2.2 mmD x 41 mmL	Drill 2.8 mmD x 29 mmL	Drill 2.8 mmD x 33 mmL	Drill 2.8 mmD x 41 mmL	Drill 3.5 mmD x 29 mmL	Drill 3.5 mmD x 33 mmL	Drill 3.5 mmD x 41 mmL	Drill 4.2 mmD x 29 mmL	Drill 4.2 mmD x 33 mmL	Drill 4.2 mmD x 41 mmL	Profile Drill 3.3 mmD x 25 mmL	Profile Drill 3.3 mmD x 33 mmL	Profile Drill 4.1 mmD x 25 mmL	Profile Drill 4.1 mmD x 33 mmL	Profile Drill 4.8 mmD x 25 mmL	Profile Drill 4.8 mmD x 33 mmL	Implant Driver NC, 3.3 mmD x 28 mmL	Implant Driver NC, 3.3 mmD x 37 mmL	Implant Driver RC, 4.1 mmD x 28 mmL	Implant Driver RC, 4.1 mmD x 37 mmL
1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2
3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3
4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4
5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5	5
6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6
7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7	7
8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8	8
9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9	9
10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10
11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11	11
12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12	12
13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13
14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14	14
15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15	15
16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16	16
17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17	17
18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18	18
19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19	19
20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20

Time to Order New Drills

\*Please consult the Reusable Instrument Lifespan Manual [ZVINST0189] for additional information. This chart is intended to be used as a reference guide only. Improper use or use in very dense bone may lead to premature wear and the need to replace the drill after fewer than 20 uses. NOTE: Surgical instruments are susceptible to damage for a variety of reasons including but not limited to: prolonged use, misuse, rough or improper handling, inadequate cleaning, etc. Care must be taken to avoid compromising the intended performance of the instrument. Please refer to P-ZBDINSTRP for complete instructions on the sterilization and care of instruments.

NC = Narrow Connection, RC = Regular Connection

For more information, visit [ZimVie.com](http://ZimVie.com)  
ZimVie US Corp LLC  
4555 Riverside Drive Palm Beach Gardens, FL 33410  
Tel: +1-561-776-6700/Fax: +1-561-776-1272  
dentalCS@ZimVie.com

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. and its affiliates; and all products are manufactured by one or more of the dental subsidiaries of ZimVie US Corp LLC (Biomet 3i, LLC, etc.) and distributed by ZimVie US Corp LLC and its authorized marketing partners. EU MDR classification: Drills and Drivers (Class IIa). The intended use of Drills is to fix prosthetic restorations and auxiliary abutments. The products here shown are not eligible for reimbursement in France. For other countries, please check the local laws and regulation. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZVINST2422-GLBL REV A 02/26 ©2026 ZimVie. All rights reserved.

