




The Evolution
of Fusion



Avenue®-T Ti




Avenue®-P Ti



Avenue®-C Ti



Avenue®-L Ti



Avenue®-A Ti

Avenue® Ti

3D Printed Fusion Solutions





Introducing **Avenue® Ti** Porous Ti Interbody System

Advancing patient care with our newest 3D printed titanium interbody platform. Avenue Ti's interbody spinal cages are created with a distinctive internal porous lattice ("net") structure. This scaffold design, structured from uniformly interconnected pores ranging from 500-700µm, along with the 6-10µm micron roughened surface topography helps to foster a cellular relevant environment for adhesion and bone ingrowth.

This range of devices have been engineered for both improved on-growth and ingrowth¹, compared to PEEK and solid titanium cages, and to comprise of a modulus of elasticity close to that of bone.

Studies have shown that the lattice structure of 3D printed titanium cages, which provide a bone comparable modulus of elasticity, is designed to withstand loading and promote fusion through providing a porous framework for bony in-growth², can sustain intra-disc height, can help to reduce the occurrence of subsidence compared to solid titanium cages³ and PEEK cages^{2,4} and are designed to provide a more evenly distributed endplate pressure under static load, compared to solid titanium and PEEK spinal cages⁴.

The Evolution of Fusion

The Avenue Ti is designed to have the following surface, structural, and anatomic features:

Porosity

- 3D printed titanium interbody spinal cages
- Balance of porosity and strength
- Engineered with an internal porous lattice structure of uniformly interconnected pores ranging in size from 500-700µm

Texture

- Microporous surface roughness of 6-10µm for potential cellular adhesion

Structure

- Able to withstand loading and promote Fusion
- Internal lattice structure which is intended to provide additional surface area for better colonisation and bone formation


Anatomic Design

- Matching patients anatomy and surgeons preferences
- Microporosity and surface roughness designed to increase friction and limit micromotion for excellent stability












Portfolio Offerings

- Static, Built-in Fixation, and Expandable options
- Wide variety of Footprints, Heights, and Lordosis Angles
- Also available in 3D Printed Tantalum





Avenue* Ti – Porous Ti Interbody Platform | 3D Printed Titanium

Description		Footprint	Lordosis	Height	
	Avenue- C Ti	ACDF	14 x 12 mm, 14 x 14 mm, 16 x 4 mm, 18 x 16 mm	5°	4 - 9 mm (1 mm increments)
	Avenue- C Fix Ti	ACDF with integrated Fixation	14 x 12 mm, 14 x 14 mm, 16 x 4 mm, 18 x 16 mm	5°	4 - 9 mm (1 mm increments)
	Avenue- P Ti	PLIF	24 x 10 mm, 29 x 10 mm	0° ⁽¹⁾ , 5°, 8° ⁽²⁾ , 14° ⁽³⁾	7 - 15 mm (1 mm increments)
	Avenue- P Exp Ti	expandable PLIF	29 x 9 mm, 32 x 9 mm	from 0° to 20° (continuous adjustability)	from 9 to 14 mm (continuous adjustability)
	Avenue- T Ti	TLIF	26 x 9 mm ⁽⁴⁾ , 29 x 9 mm, 32 x 9 mm, 32 x 10 mm	0°, 5°, 8° 15° ⁽⁵⁾	7 - 15 mm (1 mm increments)
	Avenue- T Exp Ti	expandable TLIF	32 x 9 mm, 40 x 9 mm ⁽⁶⁾	from 2° to 8° or from 6° to 14° (continuous adjustability)	from 8 to 13 mm or from 9 to 13 mm (continuous adjustability)
	Avenue- L Ti	LLIF	42 x 18 mm, 48 x 18 mm, 52 x 18 mm, 58 x 18 mm	5°, 8°, 14° ⁽⁸⁾	7 mm, 9 mm, 11 mm, 13 mm ⁽⁹⁾
	Avenue- L Fix Ti	LLIF with integrated Fixation	42 x 22 mm, 48 x 22 mm, 52 x 22 mm, 58 x 22 mm 42 x 18 mm, 48 x 18 mm, 52 x 18 mm, 58 x 18 mm	5°, 8°, 14°	7 mm, 9 mm, 11 mm, 13 mm
	Avenue- L Exp Ti	expandable LLIF	48 x 18 mm, 52 x 18 mm, 48 x 22 mm, 52 x 22 mm	from 0° to 14° (continuous adjustability)	from 8 to 13 mm (continuous adjustability)
	Avenue- A Ti	ALIF	30 x 24 mm ⁽⁷⁾ , 32 x 22 mm, 36 x 24 mm, 38 x 28 mm	8°, 14° ⁽⁹⁾	8 - 15 mm (1 mm increments) ⁽⁹⁾
	Avenue- A Fix Ti	ALIF with integrated Fixation	32 X 22 mm, 38 X 28 mm	8°, 14° ⁽⁹⁾	9 - 15 mm (1 mm increments)

(1) only in 29 mm, (2) optional in 29 mm, (3)(4)(6)(7) optional, (5) 15° are optional and available for 32 x 10mm only. (8)(9) some heights and angles are optional.
NOTE: Variations of sizes may not be available in all markets

References:

1. Rao, P.J., et al., Spine interbody implants: material selection and modification, functionalization and bioactivation of surfaces to improve osseointegration. *Orthop Surg*, 2014. 6(2): p. 81-9.
2. Chan, J.L., et al., Evolution of Bioactive Implants in Lateral Interbody Fusion. *Int J Spine Surg*, 2022. 16(S1): p. S61-S68.
3. Zhu, Y., et al., Effect of Elastic Modulus on Biomechanical Properties of Lumbar Interbody Fusion Cage. *Journal of Materials Sciences and Technology*, 2009. 25(03): p. 325-328.
4. Fogel, G., et al., Choice of Spinal Interbody Fusion Cage Material and Design Influences Subsidence and Osseointegration Performance. *World Neurosurg*, 2022. 162: p. e626-e634.

For more information, visit [ZimVie.com](https://www.zimvie.com)



BIOMET 3i Dental Iberica S.L.U.
WTC Almeda Park, Ed. 4, Planta 2
C/Tirso de Molina, 40
08940 - Cornellà de Llobregat
(Barcelona) Spain



Tsunami Medical, S.r.l.
HQ: Via E. Giorgi 27 - 41124 Modena, Italy
OHQ: Via XXV Aprile 22 - 41037 Mirandola, Italy
Phone: +39 0535 38397
Fax: +39 0535 38399

 **ZimVie**
Restoring Daily Life®

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to ZimVie, Inc. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of ZimVie. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit www.zimvie.com for additional product information. Product clearance and availability may be limited to certain countries/regions. This material is intended for health care professionals only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. CE mark in brochure is not valid unless there is a CE mark on the product label. These are medical devices that comply with the current regulation. Cages / Implants (Class IIb), Instruments / Tools (Class I). **Not for distribution in France.** ZV0931 REV B 03/24 ©2024 ZimVie, Inc. All rights reserved.

